

## Order form

BodiLink<sup>®</sup> Head Support Basic

HCPCS coding: E0955 / Cushioned headrest

E1028 / GT removable hardware

GT Lateral Hardware  
Shown on Acta-Back<sup>®</sup>

Company Name/ACCT #: \_\_\_\_\_

P.O. Number: \_\_\_\_\_

Requested By: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Ship To: \_\_\_\_\_

Client Reference: \_\_\_\_\_

For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader after saving to your local drive. Then email or print and fax to Customer Support at orders.comfort@permobil.com or 1.406.522.8563.

\*DO NOT SEND PROTECTED HEALTH INFORMATION. IT IS NOT NEEDED TO MAKE THE PRODUCT YOU ARE REQUESTING.\*

\*Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process.\*

1 BodiLink<sup>®</sup> Head support pad

- You may make one selection per column.
- Choose pad size in 1.1.
- Choose a cover material in 1.2.

## 1.1 Basic pad size

HCPCS Code: E0955 MSRP \$147.00

Length (L) refers to the actual size dimension of the support from bottom to top edge.

W		
L	6W x 5L	BL-HSP3-6W5L
L	8W x 5L	BL-HSP3-8W5L
L	10W x 6L	BL-HSP3-10W6L
L	12W x 6L	BL-HSP3-12W6L
L	14W x 6L	BL-HSP3-14W6L
L	18W x 6L	BL-HSP3-18W6L



Basic head support pads can be formed to any shape using the moldable steel insert inside the pad.

## 1.2 Basic pad cover material

**COMFORT-TEK<sup>®</sup>**

For fluid protection & an easily cleaned surface.

**STRETCH-AIR<sup>®</sup>**


For comfort & heat dissipation.

	Comfort-Tek (COMFORT-TEK)
	Stretch-Air (STRETCH-AIR)

2 BodiLink<sup>®</sup> Head support hardware

## 2.1 GT hardware, removable

HCPCS code: E1028



## Range of motion

GT model	Total # of links	Total link length	Max. anterior reach	Max. posterior reach	Max. lateral reach
LAT	5	10.5"	10"	2"	7.5"
MID	3	8.5"	8"	2"	N/A

Style	MSRP
LATERAL	\$260.00
MID-LINE	\$218.00

Part number
(BL-HS-GT3RMLAT)
(BL-HS-GT3RMMID)


## 2.2 Head support mounting option



### MOUNTING BLOCK- POSTERIOR A or B

**A:** For Acta-Back®, Acta-Embrace™, Visco-Back, and Acta-Relief with Power Mount Only

**B:** For Radius Back® or OEM wood/t-nut back supports



### MOUNTING BLOCK- QUANTUM A

For Quantum Tru-Comfort back supports



### MOUNTING BLOCK- ANTERIOR A

For AGILITY and other OEM back supports




### MOUNTING BLOCK- PERMOBIL B

For Permobil Corpus 3G Back Supports (20"H)



### MOUNTING BLOCK- PERMOBIL A

For Permobil Corpus 3G Back Supports (23"- 28"H)



### BRACKET WITH MOUNTING BLOCK- ACTA-RELIEF™

18W Shown Above  
For Non-Power Mount Only.

Mounting block OR bracket with mounting block	MSRP	Part number
Posterior A Mounting Block	\$0.00	(HS-MB-POSTA) ✓
Posterior B Mounting Block	\$0.00	(HS-MB-POSTB) ✓
Anterior A Mounting Block	\$0.00	(HS-MB-ANTA) ✓
Quantum A Mounting Block	\$101.00	(HS-MB-PW1A)
Permobil A Mounting Block	\$101.00	(HS-MB-PW2A)
Permobil B Mounting Block	\$101.00	(HS-MB-PW2B)
*Acta-Relief™ 14W Bracket with Mounting Block	\$84.00	(HS-MB-AR14W)
*Acta-Relief™ 16W Bracket with Mounting Block	\$84.00	(HS-MB-AR16W)
*Acta-Relief™ 18W Bracket with Mounting Block	\$84.00	(HS-MB-AR18W)
*Acta-Relief™ 20W Bracket with Mounting Block	\$84.00	(HS-MB-AR20W)

\*For Acta Relief™, includes mounting bracket and POSTA block.  
If you already have an Acta-Relief™ head support bracket, order POSTA Block

## 3 Head support additional mounting hardware

### 3.1 Extra mounting hardware



### MOUNTING BLOCK- POSTERIOR A or B

**A:** For Acta-Back®, Acta-Embrace™, Visco-Back, and Acta-Relief with Power Mount Only

**B:** For Radius Back® or OEM wood/t-nut back supports



### MOUNTING BLOCK- PERMOBIL A

For Permobil Corpus 3G Back Supports (23"- 28"H)




### MOUNTING BLOCK- ANTERIOR A

For AGILITY and other OEM back supports




### MOUNTING BLOCK- PERMOBIL B

For Permobil Corpus 3G Back Supports (20"H)



### MOUNTING BLOCK- QUANTUM A

For Quantum Tru-Comfort back supports



### BRACKET- ACTA-RELIEF™

18W Shown Above

Mounting block OR bracket	MSRP	Part number
Posterior A Mounting Block	\$60.00	(BL-HS-MBPOSTA)
Posterior B Mounting Block	\$60.00	(BL-HS-MBPOSTB)
Anterior A Mounting Block	\$60.00	(BL-HS-MBANTA)
Quantum A Mounting Block	\$103.00	(BL-HS-MBPW1A)
Permobil A Mounting Block	\$103.00	(BL-HS-MBPW2A)
Permobil B Mounting Block	\$103.00	(BL-HS-MBPW2B)
*Acta-Relief™ 14W Bracket	\$173.00	(BL-HS-AR14W)
*Acta-Relief™ 16W Bracket	\$173.00	(BL-HS-AR16W)
*Acta-Relief™ 18W Bracket	\$173.00	(BL-HS-AR18W)
*Acta-Relief™ 20W Bracket	\$173.00	(BL-HS-AR20W)

\*Bracket only, does not include mounting block.



Additional	MSRP	Extra part
1/4" SPACER PLATE	\$43.00/each	(BL-HS-SPCRA)