

Feather



WHO WE ARE

After 20 years of working with wheelchairs, Feather realized that some of the biggest difficulties people face when using wheelchairs is they are often heavy, expensive, and hard to maneuver. This led to the design of the Feather Chair. Allowing any user or assistant to easily operate and lift the chair when transportation or storage is needed. Having a light and collapsible chair can give people the ability to travel more and experience new places.



WHAT WE DO

Feather Mobility designs the lightest wheelchairs and scooters on the planet! Super light mobility products allow seniors and their caregivers to lift their wheelchairs in and out of a car.

Contact your McKesson Account Representative for more information

Feather Power Chair

- Just under 33 pounds
- Lightweight, foldable, and portable
- Airline-approved removable battery



Right Hand Joystick

- McKesson Number: 1224575
- MFR Number: FCP18-BK

Left Hand Joystick

- McKesson Number: 1224576
- MFR Number: FCP18-BK-L

Feather Chair

- Only 13 lbs with wheels removed
- Folds for easy storage & travel
- Simply lift into car or plane
- Supports users up to 250 lbs



- McKesson Number: 1224572
- MFR Number: EB-FCM18-BK-BKC

Feather Chair HD

- Extra-wide seat measures 22 inches
- Supports users up to 375 lbs
- Lightweight, 22 lbs wheelchair



- McKesson Number: 1224573
- MFR Number: FCM22-BK-BKC

Feather Transport with Handbrakes

- Lightweight transport chair weighs 13 pounds
- Great for travel! Easily lifts in and out of a car



- McKesson Number: 1224571
- MFR Number: EB-FCTM18-BK-BKC-HB

Feather Transport Chair

- Lightweight transport chair weighs 13 pounds
- Great for travel! Easily lifts in and out of a car



- McKesson Number: 1224570
- MFR Number: EB-FCTM18-BK-BKC

Feather Scooter

- 37-pound electric scooter
- Airline-approved with removable battery
- One-fold technology for easy travel and transport



- McKesson Number: 1224574
- MFR Number: FCS-BK