SELECTIS

Comfort Select Pressure Redistribution Bariatric Foam Mattress

The three-layered, high density foam mattress provides both comfort and optimal pressure redistribution. It is designed stronger and sturdier to accomodate the higher weight capacity.

RAISED PERIMETER

SPECS

Cover Material	Stretch, Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	3
Dedicated Heels Section	Yes
Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	600/650/700 lbs.
Warranty	Core: 1 Year Non- Prorated Cover: 1 Year Non- Prorated

PRODUCT INFO

Item #	Raised Side Rails	Product Description	Product Dimension	UOM
61056	No	Multi-Layered/ Multi Zoned Foam Mattress	42" (W) x 80" (L) x 6" (H)	1/ea
61071	Yes		42" (W) x 80" (L) x 6" (H)	1/ea
61067	No		42" (W) x 84" (L) x 6" (H)	1/ea
61069	No		48" (W) x 80" (L) x 6" (H)	1/ea
61068	No		48" (W) x 84" (L) x 6" (H)	1/ea
61070	No		54" (W) x 80" (L) x 6" (H)	1/ea

FEATURES

- Three-layered, high-density foam provides comfort and optimal pressure redistribution
- Fluid resistant, easy to clean, vapor permeable stretch mattress cover ensures maximum user comfort
- Ultra soft memory foam heel section designed to provide extremely low skin interface pressure
- Top foam layer is high resilient and die-cut providing comfort and proper pressure redistribution
- High- density foam perimeter provides support and safe patient transfer
- Meets fire code 1633