


## **Comfort Select**

Pressure Redistribution

Foam Mattress

The three-layered, high density foam mattress provides both comfort and optimal pressure redistribution. The stretch mattress cover is fluid resistant, easy to clean and vapor permeable.



Cover Material	Stretch, Fluid-Resistant Nylon Top   Waterproof Vinyl Bottom	
Mattress Layers of Foam	3	
Side Rails	Optional	
Dedicated Heels Section	Yes	
Fire Rating	CFR 16 Part 1633 Federal Fire Code	
Weight Capacity	350 lbs.	
Warranty	Core: 7 Year Non-Prorated Cover: 2 Year Non-Prorated	

## FEATURES

- Three-layered, high-density foam provides comfort and optimal pressure redistribution
- Fluid resistant, easy to clean, vapor permeable stretch mattress cover ensures maximum user comfort
- Ultra soft memory foam heel section designed to provide extremely low skin interface pressure
- Top foam layer is high resilient and die-cut providing comfort and proper pressure redistribution
- High- density foam perimeter provides support and safe patient transfer
- Meets fire code 1633

## PRODUCT INFO

Item #	Raised Side Rails	Product Description	Product Dimension	UOM
61063	No	Multi-Layered/Multi Zoned Foam Mattress	36" (W) x 76" (L) x 6" (H)	1/ea
61072	Yes	Multi-Layered/Multi Zoned Foam Mattress   Raised Perimeter and Cut-out	36" (W) x 76" (L) x 6" (H)	1/ea
61055	No	Multi-Layered/Multi Zoned Foam Mattress	36" (W) x 80" (L) x 6" (H)	1/ea
61065	Yes	Multi-Layered/Multi Zoned Foam Mattress   Raised Perimeter and Cut-out	36" (W) x 80" (L) x 6" (H)	1/ea
61064	No	Multi-Layered/Multi Zoned Foam Mattress	36" (W) x 84" (L) x 6" (H)	1/ea
61066	Yes	Multi-Layered/Multi Zoned Foam Mattress   Raised Perimeter and Cut-out	36" (W) x 84" (L) x 6" (H)	1/ea